

# Method Development of an Analytical Procedure for the Determination of Clofarabine in Pharmaceutical Formulations

Venkata Kishore<sup>1\*</sup>, V.Tejeswara Rao<sup>2</sup>, D.Aditya Deepthi<sup>3</sup>, N. Annapurna<sup>1</sup>

<sup>1</sup>Department of Chemistry, AUEC (A), Visakhapatnam.

<sup>2</sup>Department of Chemistry, MVR College, Visakhapatnam, Andhra Pradesh 530026, India.

<sup>3</sup>Department of Chemistry, St. JOSEPHS COLLEGE FOR WOMEN (A), Gnanapuram.

\*Corresponding Author: Venkata Kishore, Department of Chemistry, AUEC (A), Visakhapatnam.

Received date: November 25, 2021; Accepted date: December 08, 2021; Published date: December 13, 2021

**Citation:** V Kishore, T Rao, A Deepthi, N Annapurna. (2021). Method Development of an Analytical Procedure for the Determination of Clofarabine in Pharmaceutical Formulations. *J. Pharmaceutics and Pharmacology Research*. 4(5); DOI: [10.31579/2693-7247/060](https://doi.org/10.31579/2693-7247/060)

**Copyright:** © 2021 Venkata Kishore, This is an open access article distributed under the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

## Abstract

A novel, simple and economic high performance liquid chromatography (HPLC) method has been developed for the estimation of Clofarabine in bulk and tablet dosage form with greater precision and accuracy. The method was validated as per ICH guidelines. Validation studies demonstrated that the proposed HPLC method is simple, specific, rapid, reliable and reproducible. Hence the proposed method can be applied for the routine quality control analysis of Clofarabine in bulk and tablet dosage forms. All the components of the system are controlled using SCL-10Avp System Controller. Data acquisition was done using LC Solutions software.

**Key words:** clofarabine, RP-HPLC, method development, validation, ICH guidelines

## Introduction

For many years, the Southern Research Institute has had a programme, supported by the US National Cancer Institute, searching for new nucleoside anticancer drugs. In the early 1980s, two adenine-containing nucleosides, now known as fludarabine (Fludara; Berlex Oncology) and cladribine (Leustatin; Ortho Biotech) were in clinical trials. At the time, it was not clear whether either drug would gain approval by the FDA because some concerns were raised during preclinical and clinical development of these agents. Both drugs were susceptible to glycosidic bond cleavage with fludarabine subject to some phosphorylase cleavage and cladribine subject to both hydrolytic and enzymatic cleavage [1].

Clofarabine administered intraperitoneally had significant activity against a wide variety of human tumour xenografts implanted subcutaneously in athymic nude or severe combined immune deficiency mice [2]. Moderate to excellent sensitivity to tumour growth delays were seen in all eight human colon tumours, three out of four human renal tumours, all four non-small-cell lung tumours, and all three prostate tumours. This spectrum of widespread anticancer activity has been confirmed by other investigators in human tumour xenograft models in mice [3]. The anticancer activity of clofarabine was dose- and schedule-dependent, and greater antitumour activity was associated with more frequent administration [4]. Clofarabine is a second generation purine nucleoside analog with antineoplastic activity. Clofarabine is phosphorylated intracellularly to the cytotoxic active 5'-triphosphate metabolite, which

inhibits the enzymatic activities of ribonucleotide reductase and DNA polymerase, resulting in inhibition of DNA repair and synthesis of DNA and RNA [5-7].

Acute leukaemia is the most common paediatric cancer, with acute lymphoblastic leukaemia (ALL) and acute myelogenous leukemia (AML) being the two most common types. In the United States alone, ~2,000 children are diagnosed each year with ALL and 500 with AML [8]. Successful treatment of paediatric ALL and AML involves intensive, multi-cyclic therapy with multiple drugs that have various mechanisms of action and dosing regimens [9-10]. Such intense, cyclic treatment regimens with many different agents have reported a projected 5-year disease-free survival of 70% for paediatric patients with ALL and a complete response (CR) rate of 90% in certain forms of childhood acute leukaemia [11].

In this regard and view of the need for a suitable analytical HPLC method for routine analysis of Clofarabine in formulations. Attempts were made to develop simple, precise and accurate analytical methods for estimation of Clofarabine and extend it for their determination in formulation.

## Aim and Objective

The aim of the method is to develop an analytical procedure for the determination of Clofarabine in Pharmaceutical Formulations. The analytical procedure for determination of Assay in finished product of Clofarabine Injection, 1mg/mL is an In-House procedure.

The method shall be validated for the following parameters:

- A) **Error! Reference source not found.**
- **Error! Reference source not found.**
  - **Error! Reference source not found.**
- B) **Error! Reference source not found.**
- Interference
  - **Error! Reference source not found.**
  - Linearity
- C) **Error! Reference source not found.**
- D) **Error! Reference source not found.**
- E) **Error! Reference source not found.**
- F) Stability of Analyte in solution
- G) Filter compatibility.

System Suitability **Error! Reference source not found.**

## Method of Preparation

### Instrumentation, Chromatographic Conditions & Method:

The Chromatographic system consisted of a Shimadzu Class VP Binary pump LC-10ATvp, SIL-10ADvp Auto sampler, CTO-10Avp Column Temperature Oven, SPD-10Avp UV-Visible Detector. All the components of the system are controlled using SCL-10Avp System Controller. Data acquisition was done using LC Solutions software.

The mobile phase consisted of 85:15 (v/v) of buffer solution and acetonitrile operated on isocratic mode. The flow rate is 1.0 ml/min. Chromatographic Estimation of Clofarabine was performed on Inertsil ODS-2 (150 x 4.6) mm, 5µm column. The wavelength of detection is 263 nm. The injection volume is 25µL.

### Chromatographic conditions

A High Performance liquid chromatography equipped with UV detector and an auto sampler or its equivalent

|                      |   |                                    |
|----------------------|---|------------------------------------|
| Column               | : | Inertsil ODS-2 (150 x 4.6) mm, 5µm |
| Detection wavelength | : | 263 nm                             |
| Flow rate            | : | 1.0 mL / min                       |
| Injection volume     | : | 25µL                               |
| Run time             | : | 10 min                             |
| Column temperature   | : | 40°C                               |
| Sample cooling rack  | : | 25°C                               |

### Standards

Clofarabine working standard.

### Preparation of Buffer Solution

Dilute 1.0 mL of Glacial acetic acid in 1000 mL of water and mix.

### Preparation of Mobile phase

Mix buffer solution and acetonitrile in 85:15 (v/v) portion, sonicate well.

### Preparation of Diluents

Methanol 0.9% sodium chloride

### Preparation of 0.9% Sodium Chloride Solution

Weigh and transfer accurately 0.9 g of sodium chloride into 100 mL volumetric flask dissolve and make up to the volume with water.

### Preparation of Blank

Dilute 1.0 mL of methanol in 10 mL volumetric flask and make up to the mark with 0.9 % sodium chloride solution.

### Preparation of Standard Solution

Weigh and transfer about 30 mg of Clofarabine standard in a 50 mL volumetric flask. Add about 10 mL of methanol and sonicate to dissolve. Make up to the volume with methanol and mix well. Transfer 5 mL of above solution in to 50 mL volumetric flask, dilute and make up to the volume with 0.9% sodium chloride solution.

**Note:** Standard solutions are stable up to 48 hrs at both room temperature and 2-8°C.

### Preparation of Test Solution

Transfer 3.0 mL of sample solution into 50 mL volumetric flask. Add 5.0 mL methanol and 30 ml of 0.9% sodium chloride and mix well. Dilute up to the mark with 0.9% sodium chloride solution.

**Note:** Test solutions are stable up to 48 hrs at both room temperature and 2-8°C.

### Procedure

Separately inject each solution into the chromatographic system in the following order.

|                              |   |                  |
|------------------------------|---|------------------|
| Blank injection              | - | Single           |
| Standard solution injections | - | Five             |
| Test solution                | - | Two injections   |
| Standard solution bracketing | - | Single injection |

## Results and Discussion

### Experimental:

#### **Error! Reference source not found.**

Intermediate precision expresses within-laboratories variations such as different days, different analysts, different columns, different equipment's etc. Ruggedness incorporates the concept described under the terms "Intermediate Precision" as defined in USP <1225>.

Intermediate precision is established by doing same exercise as system and method precision by different analyst on different day using different column and different equipment. The same Lot/Batch of standard and sample were used within the laboratory. The results of Intermediate Precision are tabulated in below table. Comparison of Method Precision and Intermediate Precision result is summarized in the table below.

### Results of Intermediate Precision (Sample Solution)

| Sample #     | % Assay |
|--------------|---------|
| 1            | 101.1   |
| 2            | 101.1   |
| 3            | 101.1   |
| 4            | 101.1   |
| 5            | 100.2   |
| 6            | 100.2   |
| <b>Mean</b>  | 101.1   |
| <b>% RSD</b> | 0.05    |

### Comparison of Method Precision and Intermediate Precision Results

| Parameter   | Method Precision       | Intermediate Precision |
|---|------------------------|------------------------|
| Analyst   | Analyst 1              | Analyst 2              |
| HPLC ID.  | EP-QCI-012             | EP-QCI-013             |
| Column ID.  | HPLCC-008              | HPLCC-009              |
| Comparison of Method Precision and Intermediate precision |                        |                        |
| Sample #  | % Assay of Clofarabine |                        |
|   | Method Precision       | Intermediate Precision |
| 1   | 101.9                  | 101.1                  |
| 2   | 101.9                  | 101.1                  |
| 3   | 101.9                  | 101.1                  |
| 4   | 101.9                  | 101.1                  |
| 5   | 101.9                  | 101.2                  |
| 6   | 102.1                  | 101.2                  |
| <b>Mean</b>   | 101.9                  | 101.1                  |
| <b>% RSD</b>  | 0.02                   | 0.05                   |
| <b>Overall Mean (12 samples)</b>                          | 101.5                  |                        |
| <b>Cumulative % RSD (12 samples)</b>                      | 0.42                   |                        |

#### Acceptance Criteria:

1. The relative standard deviation of results obtained from six sample preparations should not be more than 2.0%
2. The cumulative relative standard deviation of method precision and intermediate precision results obtained from twelve sample (6 methods precision and 6 intermediate precision) preparations should not be more than 2.0%.

#### Conclusion:

The result meets the acceptance criteria and found comparable, indicates that the method is precise and rugged with respect to analyst to analyst, day to day, column to column and equipment to equipment for its intended use.

Error! Reference source not found.

The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage. Robustness study is performed by analyzing the standard at different conditions. The results obtained with altered conditions are compared against results obtained under normal chromatographic conditions.


#### Variation in Flow Rate ( $\pm 0.2$ mL/min.)

The standard was carried out by varying the flow rate of mobile phase to 0.8 mL/min. and 1.2 mL/min. in place of actual flow rate 1.0 mL/min. The results are summarized in the below table. Typical chromatogram of Robustness for variation in flow rate (0.8 mL/min and 1.2 mL/min) is exhibited below.

#### Results of robustness -Variation in flow rate for Clofarabine

| Injection #               | Flow Rate<br>0.8 mL/min. | Actual Flow Rate<br>1.0 mL/min. | Flow Rate<br>1.2 mL/min. |
|---------------------------|--------------------------|---------------------------------|--------------------------|
| 1                         | 5449522                  | 4365531                         | 3613880                  |
| 2                         | 5448777                  | 4366135                         | 3612625                  |
| 3                         | 5451470                  | 4366853                         | 3613455                  |
| 4                         | 5444058                  | 4364875                         | 3613092                  |
| 5                         | 5449326                  | 4366851                         | 3614659                  |
| <b>Mean</b>               | 5448631                  | 4366049                         | 3613542                  |
| <b>% RSD</b>              | 0.05                     | 0.02                            | 0.02                     |
| <b>Tailing factor</b>     | 1.1                      | 1.1                             | 1.0                      |
| <b>Theoretical plates</b> | 6243                     | 4943                            | 3940                     |


**Chromatogram of Robustness for variation in flow rate (0.8 mL/min)**



**Peak Results**

| Name          | RT    | Area    | % Area | Height (μV) | USP Plate Count | USP Tailing | Int Type |
|---------------|-------|---------|--------|-------------|-----------------|-------------|----------|
| 1 Clofarabine | 4.536 | 5449522 | 100.00 | 633436      | 6296            | 1.1         | BB       |

**Chromatogram of Robustness for variation in flow rate (1.2 mL/min)**



**Peak Results**

| Name          | RT    | Area    | % Area | Height (μV) | USP Plate Count | USP Tailing | Int Type |
|---------------|-------|---------|--------|-------------|-----------------|-------------|----------|
| 1 Clofarabine | 3.083 | 3613880 | 100.00 | 497328      | 3911            | 1.0         | BB       |

**Variation in Column Oven Temperature (± 2°C)**


The standard was carried out by varying the column oven temperature of 38°C and 42°C in place of actual column oven temperature 40°C. The results are summarized in the below table. Chromatogram of Robustness

for variation in Column Oven Temperature (38°C and 42°C) is exhibited below.

**Results of Robustness-Variation in Column oven Temperature for Clofarabine**

| Injection #               | Column Oven Temperature 38°C | Actual Column Oven Temperature 40°C | Column Oven Temperature 42°C |
|---------------------------|------------------------------|-------------------------------------|------------------------------|
| 1                         | 4356727                      | 4365531                             | 4354231                      |
| 2                         | 4357687                      | 4366135                             | 4357138                      |
| 3                         | 4359938                      | 4366853                             | 4354543                      |
| 4                         | 4357320                      | 4364875                             | 4352076                      |
| 5                         | 4358888                      | 4366851                             | 4361549                      |
| <b>Mean</b>               | 4358112                      | 4366049                             | 4401322                      |
| <b>% RSD</b>              | 0.03                         | 0.02                                | 0.08                         |
| <b>Tailing factor</b>     | 1.1                          | 1.1                                 | 1.1                          |
| <b>Theoretical plates</b> | 4914                         | 4943                                | 4875                         |


**Chromatogram of Robustness for variation in Column Oven Temperature (38°C)**



**Peak Results**

| Name          | RT    | Area    | % Area | Height (µV) | USP Plate Count | USP Tailing | Int Type |
|---------------|-------|---------|--------|-------------|-----------------|-------------|----------|
| 1 Clofarabine | 3.720 | 4356727 | 100.00 | 549731      | 4911            | 1.1         | BB       |

**Chromatogram of Robustness for variation in Column Oven Temperature (42°C)**



**Peak Results**

| Name          | RT    | Area    | % Area | Height (µV) | USP Plate Count | USP Tailing | Int Type |
|---------------|-------|---------|--------|-------------|-----------------|-------------|----------|
| 1 Clofarabine | 3.612 | 4354231 | 100.00 | 566163      | 4948            | 1.1         | BB       |

**Variation in Organic composition (142.5 and 157.5)**


The standard was carried out by varying the Organic composition 142.5 mL and 157.5 mL in place of actual the 150mL. The results are summarized in the below table. Chromatogram of Robustness for

variation in the Organic composition (142.5 mL and 157.5 mL) is exhibited below.

**Results of Robustness-Variation in pH of Buffer for Clofarabine**

| Injection #               | Low Organic composition 142.5 mL | Actual Organic composition 150 mL | High Organic composition 157.5 mL |
|---------------------------|----------------------------------|-----------------------------------|-----------------------------------|
| 1                         | 4360826                          | 4365531                           | 4363298                           |
| 2                         | 4359986                          | 4366135                           | 4360568                           |
| 3                         | 4362608                          | 4366853                           | 4359309                           |
| 4                         | 4360929                          | 4364875                           | 4360570                           |
| 5                         | 4361818                          | 4366851                           | 4358112                           |
| <b>Mean</b>               | 4361234                          | 4366049                           | 4360371                           |
| <b>% RSD</b>              | 0.02                             | 0.02                              | 0.04                              |
| <b>Tailing factor</b>     | 1.1                              | 1.1                               | 1.1                               |
| <b>Theoretical plates</b> | 5489                             | 4943                              | 4545                              |


Chromatogram of Robustness for variation in organic composition of mobile phase (Low organic).



Peak Results

| Name          | RT    | Area    | % Area | Height (μV) | USP Plate Count | USP Tailing | Int Type |
|---------------|-------|---------|--------|-------------|-----------------|-------------|----------|
| 1 Clofarabine | 3.990 | 4360826 | 100.00 | 539588      | 5489            | 1.1         | BB       |

Chromatogram of Robustness for variation in organic composition of mobile phase (High organic).



Peak Results

| Name          | RT    | Area    | % Area | Height (μV) | USP Plate Count | USP Tailing | Int Type |
|---------------|-------|---------|--------|-------------|-----------------|-------------|----------|
| 1 Clofarabine | 3.459 | 4363298 | 100.00 | 569667      | 4611            | 1.1         | BB       |

**Acceptance Criteria:**

The System suitability defined in test procedure should meet in each condition.

- The Tailing factor for Clofarabine should be NMT 2.0.
- The relative standard deviation for Clofarabine peak from five replicate injections of standard solution should be NMT 2.0 %.

- The theoretical plates for Clofarabine peak in standard solution should be not less than 3000.

**System suitability of overall validation study**

The System suitability is an integral part of analytical procedure. The tests are based on the concept that the equipment, analytical operations and samples to be analyzed constitute an integral system that can be evaluated as such. The system suitability results are tabulated in the below Table.

**System Suitability of Overall Validation Study**

| Parameter  |      | Tailing Factor | Theoretical plates |
|--|------|----------------|--------------------|
| System Suitability/ System Precision                 | 0.02 | 1.1            | 4943               |
| Specificity by diluent, placebo and known impurities | 0.02 | 1.1            | 4943               |
| Specificity by Forced degradation                    | 0.02 | 1.1            | 4983               |
| Specificity by Forced degradation ( Alkali)          | 0.03 | 1.1            | 5047               |
| Specificity by Forced degradation ( UV)              | 0.02 | 1.1            | 4983               |
| Linearity  | 0.02 | 1.1            | 4943               |
| Method Precision                                     | 0.02 | 1.1            | 4943               |
| Intermediate Precision                               | 0.04 | 1.1            | 6303               |
| Accuracy (Recovery)                                  | 0.02 | 1.1            | 4943               |
| Robustness-Flow rate: 0.8mL/minute                   | 0.05 | 1.1            | 6243               |
| Robustness-Flow rate: 1.2mL/minute                   | 0.02 | 1.0            | 3940               |
| Robustness-Column oven temperature: 38°C             | 0.03 | 1.1            | 4914               |
| Robustness-Column oven temperature: 42°C             | 0.08 | 1.1            | 4875               |
| Robustness-Low organic composition(142.5 mL)         | 0.02 | 1.1            | 5489               |
| Robustness-High organic composition(157.5 mL)        | 0.04 | 1.1            | 4545               |
| Stability of Analyte in Solution (Initial)           | 0.02 | 1.1            | 4943               |
| Stability of Analyte in Solution (24 Hours)          | 0.02 | 1.1            | 4943               |
| Stability of Analyte in Solution (48 Hours)          | 0.02 | 1.1            | 5489               |
| <b>Minimum</b>                                       | 0.02 | 1.0            | 3940               |
| <b>Maximum</b>                                       | 0.08 | 1.1            | 6243               |
| <b>Average</b>                                       | 0.03 | 1.1            | 5078               |

**Conclusion**

The analytical procedure for Assay is validated and found suitable for its intended use and it meets the acceptance criteria for:

**Specificity:**

No Interference should be observed at the retention time of peak in the chromatograms obtained from the diluent and the placebo solution.

**Forced Degradation:**

The method is specific and stability indicating for its intended use.

**Linearity:**

The analytical procedure is linear within the concentration range from 50 % to 150 % (30.24µg/mL to 90.72µg/mL) for Clofarabine peak.

**Intermediate Precision:**

The method is precise and rugged with respect to analyst to analyst, day to day, column to column and equipment to equipment for its intended use.

**Accuracy:**

The analytical test procedure is accurate for its intended use.

**Robustness:**

The test method is robust enough as demonstrated by altering the Flow rate, Column oven temperature and Organic composition.

**Stability of analyte in solution:**

The Standard solution is stable up to 48 hours and sample solution is stable up to 48 hours at room temperature and 2-8°C for Clofarabine peak.

The data for each validation characteristic described in this report meets the acceptance criteria with respect to Specificity, Forced degradation, Stability of analyte in solution, Linearity, Precision, Intermediate Precision, Accuracy and Robustness.

The validation results reveal that the analytical procedure is suitable for determination of Assay in ClofarabineInjection, 1mg/mL. The method is stability indicating for determination of Assay of Clofarabine in ClofarabineInjection, 1mg/mL.

## References

1. Lindemalm S, Liliemark J, Juliusson J, Larsson R. & Albertioni F. (2004). Cytotoxicity and pharmacokinetics of cladribine metabolite, 2-chloroadenine, in patients with leukemia. *Cancer Lett.* 210:171-177.
2. Waud, W. R., Schmid, S. M., Montgomery, J. A. & Secrist, J. A. (2000). Antitumor activity of 2-chloro-9-(2-deoxy-2-fluoro-β-D-arabinofuranosyl) adenine (Cl-F-ara-A). *Nucleosides Nucleotides Nucleic Acids.* 19:447-460.
3. Takahashi T, Kanazawa J, Akinaga S, Tamaoki T & Okabe M. (1999). Antitumor activity of 2-chloro-9-(2-deoxy-2-fluoro-β-D-arabinofuranosyl) adenine, a novel deoxyadenosine analog, against human colon tumor xenografts by oral administration. *Cancer Chemother. Pharmacol.* 43:233-240.
4. Stephenson, K. et al. (2003). Correlation between frequency of administration and efficacy of clofarabine in the H460 human non-small cell lung tumor xenograft model. *Proc. Am. Assoc. Cancer Res.* 44.
5. DeGennaro LJ, Raetz E. (2014). Acute Lymphoblastic Leukemia.
6. Venkata NR, Jalandhar D, Gnanadey G, Bandari R, Manoi P. (2013). Development of supercritical Fluid (Carbon dioxide) based ultra-performance convergence chromatographic stability indicating assay method for determination of clofarabine in injection. *Analytical methods.* 5:7008-7013.
7. Takahia Y, Yamauchi T, Rie N, Takanori U (2011). Determination of Clofarabine triphosphate concentration in leukemia cell using <sup>32</sup>P. *Anal. Biochem.* 414:186-190.
8. Pui C. H. (1995). Childhood leukemias. *N. Engl. J. Med.* 332:1618-1630.
9. Nachman J. B. et al. (1998). Augmented post-induction therapy for children with high-risk acute lymphoblastic leukemia and a slow response to initial therapy. *N. Engl. J. Med.* 338, 1663-1671.
10. Nachman J. et al. (1997). Augmented Berlin–Frankfurt–Munster therapy abrogates the adverse prognostic significance of slow early response to induction chemotherapy for children and adolescents with acute lymphoblastic leukemia and unfavorable presenting features: a report from the Children’s Cancer Group. *J. Clin. Oncol.* 15:2222-2230.
11. Kersey J. H. (1997). Fifty years of studies of the biology and therapy of childhood leukemia. *Blood.* 90:4243-4251.



This work is licensed under Creative Commons Attribution 4.0 License

To Submit Your Article Click Here:

[Submit Manuscript](#)

DOI: [10.31579/2693-7247/060](https://doi.org/10.31579/2693-7247/060)

### Ready to submit your research? Choose Auctores and benefit from:

- fast, convenient online submission
- rigorous peer review by experienced research in your field
- rapid publication on acceptance
- authors retain copyrights
- unique DOI for all articles
- immediate, unrestricted online access

At Auctores, research is always in progress.

Learn more <https://auctoresonline.org/journals/pharmaceutics-and-pharmacology-research>